

C16orf80-chromosome 16 open reading frame 80 Gene

Reference: PTXBC005152

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0468 protein C16orf80; C16orf80; EVORF; GTL3; fSAP23; cilia- and flagella-associated protein 20; basal body up-regulated protein 22; flagellar associated protein 20 homolog; functional spliceosome-associated protein 23; gene trap locus 3; transcription		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 80
GENE ID	29105
ORIGIN SPECIES	Human
ACCESSIONS	BC005152
VECTOR	pDONR223
NCBI SYMBOL	C16orf80
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcaaaaacacgttcagagcggcttcctccatctctacagcatcgagcagcctctgcaaatctgggacaaaaaggtacggaatggccacatcaaaagaatcactgataatgacatccagtccttggtgctagagattgaaggacaaatgtaagcaccacatatcacatg
ccctgcagacccaagaagacgtgggaattaaactcttctctgtcatgattatcaaaaacctgaagaagtattttaccttcgaagtcagggtactagatgacaagaatgtcgtcgtcgttctgggcaagtaactaccagagcaccacccgggtcaaaccttcatctgcacatgccca
tgcggctggatgacggctggaaccagattcagttcaactgctagacttcacacggcgagcatagcggcacaattatcatcgagaccctcagagtcagatccatgcaaatgtcgcacccgacgggtttacttctcagacagacttactcagaagatgagctccggcagagttcaactgt
atctccagttcagaacaggcaagcaataa