



Reference: PTXBC005158

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	endosomal/lysosomal potassium channel TMEM175; hTMEM175; transmembrane protein 175		

For research use only.

GENE DESCRIPTION	transmembrane protein 175
GENE ID	84286
ORIGIN SPECIES	Human
ACCESSIONS	BC005158
VECTOR	pDONR223
NCBI SYMBOL	TMEM175
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]