



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

MRPS18C-mitochondrial ribosomal protein S18C Gene

Reference: PTXBC005186

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-134; MRP-S18-1; MRP-S18-c; MRPS18-1; S18mt-c; mrps18-c; 28S ribosomal protein S18c, mitochondrial; 28S ribosomal protein S18-1, mitochondrial; mitochondrial ribosomal protein S18-1; mitochondrial ribosomal protein S18C		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S18C
GENE ID	51023
ORIGIN SPECIES	Human
ACCESSIONS	BC005186
VECTOR	pDONR223
NCBI SYMBOL	MRPS18C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccctgctgctgttgcgtttccgggtgctataggaggagaagaagtgtacacacttggtaacggctctgctcagccttacacatcccccggactcacacggctcttggagaagagggtgttcacacaggtattccagcaatgaggactgccatftcaatggaaaactctataaagaacctt
taagaagaatgtctgtgtggaagacatgtatgaataagaattgacatcttctgccacgtttgttctccactctgcatgacalltattggaaggacattacagctctttgtggaagaacacagaagaatacacaaaagcaattagagagctcaaatatggggttatgccaggtatcacaca
agcattctcattctacacgaagcctataaattgttaatacagatatacggaaata

More info: Please contact us for any question: contact@proteogenix.fr