

MED4-mediator complex subunit 4 Gene

Reference: PTXBC005189

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARC36; DRIP36; HSPC126; TRAP36; VDRIP; mediator of RNA polymerase II transcription subunit 4; TRAP/SMCC/PC2 subunit p36; activator-recruited cofactor 36 kDa component; mediator, 34-kD subunit, homolog; vitamin D receptor-interacting protein, 36-kD; vitami		

For research use only.

GENE DESCRIPTION	mediator complex subunit 4
GENE ID	29079
ORIGIN SPECIES	Human
ACCESSIONS	BC005189
VECTOR	pDONR223
NCBI SYMBOL	MED4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcgtcttcgagtgagaggaaggagcggctgggagggcgggttgggagtgccgggtggttaacagcacacgagagcggctgctgtcgcgttgaggacttgagggtcctgtctaggggaacttatagaatgctggcaatttcaagaaccagaagttgttacaggctgg
agaggaaccaggtcctggagttgttaattcaccgagatgggaatttcaagaactaatgaaattggcacttaacagggaaaaaattcatgaaatgcaagtttagaaaaagaagtagagaagagagacggtgatattcagcagctacaaaaacagctaaaggagcagaacaaata
ctggcaacagctgtttaccaagcgaaggagaactcaagtcataagaaaaagcaaaaaaggtgctatcctctgaagaaataaagtagatgcacatagatcagtcgaagtaagctgtatgtgctccactgacctgggtccaggggacccccggagaccctaccaactgatttaga
gatgagaagtggttactgggtcagatgaacatcttcactaatggcgtgaatggcattaccaggagatgcactgcagcaggaagattgccagatgtctgtccacagatccatggcagtcgaatgacatgtcgaatatgttaccacaaatcatagtagtacttttggga
acctcctgggcataataaagaagatgaagatgatgtagagattatgtcaacggactcctcaagcagtagtagtgagctgattga