

FKBP14-FK506 binding protein 14, 22 kDa Gene

Reference: PTXBC005206

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP14; peptidyl-prolyl cis-trans isomerase FKBP14; EDSKMH; FKBP22; IPBP12; 22 kDa FK506-binding protein; 22 kDa FKBP; FK506 binding protein 14, 22 kDa; FKBP-22; rotamase; FK506 binding protein 14		

For research use only.

GENE DESCRIPTION	FK506 binding protein 14, 22 kDa
GENE ID	55033
ORIGIN SPECIES	Human
ACCESSIONS	BC005206
VECTOR	pDONR223
NCBI SYMBOL	FKBP14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctttctgtggaacgcggtctgactctgttcactctttgattgggcttgcctgaaccagaagtgaattgaagttctccagaagccatcatctgccatcgcaagaccaaggaggaggattgatgttggtccactatgaaggctacttagaaaaggacggctcttattt
cactccactcacaacataacaatggtcagccatttgggttacctgggcatctggaggctctcaagggtgggaccagggtgaaaggaaatgtgtgtaggagagaagagaagctcatctcctctgctctgggctatggaaaaggaaaaggtaaaattccccagaaagtac
actgatatataatgatctcctggagattcgaaatggaccaagatcccatgaatcattccaagaatggatcttaatgatgactggaactctctaaagatgagggttaagcatalttaagaaggagtttgaaaaacatggtgctggtggtgaatgaaatcatcatgatgcttgggtggaggat
attttgataaagaagatgaagacaagatgggttatatctgccagagaattacatataaacacgatgagttatag