

STK4-serine/threonine kinase 4 Gene

Reference: PTXBC005231

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KRS2; MST1; YSK3; serine/threonine-protein kinase 4; STE20-like kinase MST1; kinase responsive to stress 2; mammalian STE20-like protein kinase 1; mammalian sterile 20-like 1; serine/threonine-protein kinase Krs-2; serine/threonine kinase 4		

For research use only.

GENE DESCRIPTION	serine/threonine kinase 4
GENE ID	6789
ORIGIN SPECIES	Human
ACCESSIONS	BC005231
VECTOR	pDONR223
NCBI SYMBOL	STK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacgggtacagctgaggaaccgcccgcgcgcagctgaaaaagttgatgaagatagtttaaccaacaaccagaagaagtatttgatgtcttagagaacttgagaagggtga