

C1D-nuclear DNA-binding protein Gene

Reference: PTXBC005235

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1D nuclear receptor corepressor; C1D nuclear receptor co-repressor; C1D DNA-binding protein; nuclear nucleic acid-binding protein C1D; LRP1; Rrp47; SUN-CoR; SUNCOR; hC1D; nuclear DNA-binding protein; small unique nuclear receptor co-repressor; small uniq		

For research use only.

GENE DESCRIPTION	nuclear DNA-binding protein
GENE ID	10438
ORIGIN SPECIES	Human
ACCESSIONS	BC005235
VECTOR	pDONR223
NCBI SYMBOL	C1D
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggtgaagaaataatgaagactatccagtagaattcacgagatfittgcagcgtttgagaattccattggtcgtgtggatgagatgctgaagaccatgatgtctgttctagaatgagttgtgcagaagttggatccactgaacaagcaaaagtggattgtttctgcatacacatta
aattcaatgtttgggttatttggcaaccaagagagtaatcctaaaggaacatccagtaaaacaggaattgaaagaatcagagatataatgaacagagtcaaggaataacagacaagaaaaaggctggcaagctggacagaggtgcagctcaagattgtaaaaaatgccctctggaa
ccaaaatcgaaaaatgcatcaaaagtccaataaaggaaaaagcaaaagttaa