



Reference: PTXBC005237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOBK13; 2C4D; CGI-95; MOB1; MOB3; PHOCN; PREI3; MOB-like protein phocein; MOB1, Mps One Binder kinase activator-like 3; Mps One Binder kinase activator-like 3; class II mMOB1; mob1 homolog 3; phocein, Mob-like protein; preimplantation protein 3; MOB famil		

GENE DESCRIPTION	MOB1, Mps One Binder kinase activator-like 3 (yeast)
GENE ID	25843
ORIGIN SPECIES	Human
ACCESSIONS	BC005237
VECTOR	pDONR223
NCBI SYMBOL	MOBKL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgctcatgcgcggagggggcgcgcagctgctgaggcggaaagcgcgggaccacgagccaggcttctataaattggcctgatgaatcctttgatgaatggacagctacactagctgttcaacagctatattcaacagaacataagacagattgctccaattatgacaaaattcttgaaacct
gaagcccaagatgaagggtgtgtggaagtatgaacatttaaggcagtttgccttgagcttaaatgagactctgttcaaaactcagagtggaatgccatccagatacttgcatcacaatgacagcaactgaaacaattggattttctgtgtgacgtcataaaactcaaaagagtgctcctctatagact
atactagacacacactgatgtgtgacatcttctctgaatagcaataaataattttcccgacgggttagcataaaggaaatcctgttagcgaaactaggatcagtatgccgttagctattacagaatttttccatcgtctattttccatcgcgcgcagatattttgagaattatgaaatgaacaactttttg
ctatcccttttactacaattgtatgatacaatacaattgtatgctcaagagataactctgtatcaaaattttagaagaaggaaatcagacaagattctctctggcgaaagtgaagactga