

HSPB11-heat shock protein family B (small), member 11 Gene

Reference: PTXBC005245

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf41; HSPCO34; IFT25; intraflagellar transport protein 25 homolog; heat shock protein beta-11; heat shock protein family B (small), member 11; intraflagellar transport 25 homolog; placental protein 25; heat shock protein family B (small) member 11		

For research use only.

GENE DESCRIPTION	heat shock protein family B (small), member 11
GENE ID	51668
ORIGIN SPECIES	Human
ACCESSIONS	BC005245
VECTOR	pDONR223
NCBI SYMBOL	HSPB11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaaaaattgatctctgtctgagctctgaagggtccgaagtgtatttagctacatcaagtgtgaaaaacaccacacgtgaaatatcattgatgggaatccagaacgtttggaccaccacaggaatgttcccaggaattcattattgtttccacaacatgtaaggattgaaaggctgt
aatccaaagttaattgtacagacctgaagattgaaaaaacgctctaaagagccagttgattttgagcaatggattgaaaaagatttggtacacacagaggggcagcttcaaatgaagaaattgtggcacatgatggctcgcgtacttacttgagattcattattgtatcagcctttgatcatt
tgcatctgtgcatagcgtttctgcagaaggaacagtagtctcaaatctttctcataa