

SLN-sarcolipin Gene

Reference: PTXBC005261

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sarcolipin		
For research use only.			

GENE DESCRIPTION	sarcolipin
GENE ID	6588
ORIGIN SPECIES	Human
ACCESSIONS	BC005261
VECTOR	pDONR223
NCBI SYMBOL	SLN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcctgggattgactgagatgtccggagctgccgtctatgccctgagaccccaactgtctcattgtcacaggatgccattctccatccgagggcacctgtga