

RPA3-replication protein A3, 14kDa Gene

Reference: PTXBC005264

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	REPA3; RP-A p14; replication protein A 14 kDa subunit; RF-A protein 3; replication factor A protein 3; replication protein A3, 14kDa; replication protein A3		

For research use only.

GENE DESCRIPTION	replication protein A3, 14kDa
GENE ID	6119
ORIGIN SPECIES	Human
ACCESSIONS	BC005264
VECTOR	pDONR223
NCBI SYMBOL	RPA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggacatgatggacttggcaggtcgccatcaacggcgcatgctagctcaattcatcgacaagcctgtctgcttctgtagggaggtggaaaagattcatccaccggaaaaatgtttattcttcagatggagaaggaaaaatggaccatcagttgatggaacccttgatgaa
gaaatctctggaattgtggaagtgttggaaagagtaaccgccaaggccaccatctgtgtacatcttatgtccagtttaagaagatagccatctttgatcttggactttacaatgaagctgtgaaaattatccatgacttcctcagtttatcctttaggattgtgcaacatgattga