

NDUFS4-NADH dehydrogenase (ubiquinone) Fe-S protein 4, 18kDa (NADH-coenzyme Q reductase) Gene

Reference: PTXBC005270

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AQDQ; CI-18; CI-18 kDa; CI-AQDQ; NADH dehydrogenase [ubiquinone] iron-sulfur protein 4, mitochondrial; NADH dehydrogenase (ubiquinone) Fe-S protein 4, 18kDa (NADH-coenzyme Q reductase); NADH-ubiquinone oxidoreductase 18 kDa subunit; complex I 18kDa subunit		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) Fe-S protein 4, 18kDa (NADH-coenzyme Q reductase)
GENE ID	4724
ORIGIN SPECIES	Human
ACCESSIONS	BC005270
VECTOR	pDONR223
NCBI SYMBOL	NDUFS4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggtctcaatgtcagtggttactgaggcagacgttggcgggagagggcagtggtgtgagctgccctttccgtttccagggttcgaccagggtggtgaggacttctcatggagattggcacaggaccagactcaagacacacactcatacagttgatgaaaaattggatat
cactactttaactggcgttccagaagagcatataaaaactagaaaagtcaggatcttggctctgctcgaataacatgcagctctggagtaaacacacaaagaatggaagatggagtttgataccaggagcagatgggaaaaatccttgatgggttgggcatcaacggctgatcccttatcc
aacatgggttctaacttcagttactaaagaagatgcagtttcttgcagaaaaaatggatggagctatgacattgaagagagggaaggttccaaacccaagtccttatgggtgcaaacctttcttgggacaaaagaacaagagatccacaaaatag