

ZFAND6-zinc finger, AN1-type domain 6 Gene

Reference: PTXBC005283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZA20D3; ZFAND5B; AN1-type zinc finger protein 6; protein associated with PRK1; zinc finger, A20 domain containing 3; zinc finger, AN1-type domain 6; zinc finger AN1-type containing 6		

For research use only.

GENE DESCRIPTION	zinc finger, AN1-type domain 6
GENE ID	54469
ORIGIN SPECIES	Human
ACCESSIONS	BC005283
VECTOR	pDONR223
NCBI SYMBOL	ZFAND6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcaagaaataacacagccaagtgcctatgctttgtccactggctgtggattttatggaaacctctgacaaatggcatgtgttcagtagctataaagaacatcttcaagacagaatagtagtaaggtagaataagccacctgcaacctctgtagtagtctgtctgaatcttacc
agttcaatgcacagatggcagtggtccagaagcccagtcagcattagactctacatcttcatctatgcagccagccctgtatcaaatcagtcactttatcagaatctgtgacatcttctcaattggacagtagcatctgtggacaaagcagtagctgaaacagaagatgtgcaggcttcagtagtca
gacacagcacagcagccatctgaagagcaaaagcagtagtcttggaaaaaccgaaacaaaaaagaatcgctgtttcatgtgcaggagaagaagtggaacttactgggttggaatgccggtgtggaaatgtttactgtggtgtacacggtactcagatgtacacaattgctcttacaattacaag
ccgatgctgtgagaaatcagaaagaaaatccagtagttgtgtggaagaatccaaagatttga