

42628-15 kDa selenoprotein Gene

Reference: PTXBC005294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEP15; 15 kDa selenoprotein; selenoprotein F		
For research use only.			

GENE DESCRIPTION	15 kDa selenoprotein
GENE ID	9403
ORIGIN SPECIES	Human
ACCESSIONS	BC005294
VECTOR	pDONR223
NCBI SYMBOL	15-sept.
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctggccgagtggtgtctgtggtccggtttgggtctacggtgtgtgtggcgactgtgtctcaagcgggtgtctgttttggggcagagtttcatcggaggcatgcagagagtaggctttttagcaactgtcttgcagctcttgcagcttctcggacagttcaacctgtctcagct
ggatcctgattgcagagatgctgtcaggaggaagcacaatttgaacacaaaagctgtatgcaggagctattcttgaagtgtgtggaataatgggaaggttcctcaagtcgaagctttttaggagtgataaacccaaactgttcagaggactgcaaatcaagtatgccgtgttcag
accctgtattaaagcttttggacgacaatgggaacattgtgaagaactgagcattctcaaatggaacacagacagtgtagaagaattcctgagtgaaaagttggaacgcatataa