

FXYP2-FXYP domain containing ion transport regulator 2 Gene

Reference: PTXBC005302

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP1G1; HOMG2; sodium/potassium-transporting ATPase subunit gamma; ATPase, Na ⁺ /K ⁺ transporting, gamma 1 polypeptide; Na ⁽⁺⁾ /K ⁽⁺⁾ ATPase subunit gamma; sodium pump gamma chain; FXYP domain containing ion transport regulator 2		

For research use only.

GENE DESCRIPTION	FXYP domain containing ion transport regulator 2
GENE ID	486
ORIGIN SPECIES	Human
ACCESSIONS	BC005302
VECTOR	pDONR223
NCBI SYMBOL	FXYP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaggtgtgtacgtggcggcagcccaagggggacgtggaccgttctactatgactatgagaccgttcgaatggggcctgactctcgtggactggccttcacgtggggctcctcctcctcagcagaagattccgtgtgggggcaataagaagcgcaggcacaatcaatgaagatgagccgtaa