

TGFA-transforming growth factor, alpha Gene

Reference: PTXBC005308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TFGA; protransforming growth factor alpha; TGF-alpha		
For research use only.			

GENE DESCRIPTION	transforming growth factor, alpha
GENE ID	7039
ORIGIN SPECIES	Human
ACCESSIONS	BC005308
VECTOR	pDONR223
NCBI SYMBOL	TGFA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtccctcggctggacagctcgcctgttgcgtctgggtattgtgttgctgctgccccttgagaaacagcagtcctccctgagtgaccgcccgtggctgcagcagtggtgtccattttaatgactgccagattccacactcagttctgtccatggaacctgcaggtttt
gggtcaggaggaagccagcatgtgtctgcatctcgggtacgttgggtgacgctgtgagcatcgaggacctcctgcccgtggctgcccagcagaagaagcaggccatcaccgcttgggtgtgtctccatcgtggcctggtgtcttcatcatcatgtgtgtgatacactgt
gccaggtccgaaacactgtgagtggtgccggccctcatctgccggcacgagaagccagcgcctcctgaagggaagaaccgcttgcctcactcagaaacagtggtctga