



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KLK1-kallikrein 1 Gene

Reference: PTXBC005313

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KLKR; Klk6; hK1; kallikrein-1; glandular kallikrein 1; kallikrein 1, renal/pancreas/salivary; kallikrein serine protease 1; kidney/pancreas/salivary gland kallikrein; tissue kallikrein; kallikrein 1		

For research use only.

GENE DESCRIPTION	kallikrein 1
GENE ID	3816
ORIGIN SPECIES	Human
ACCESSIONS	BC005313
VECTOR	pDONR223
NCBI SYMBOL	KLK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgcttgcctgctgtgctcctgcgcctctccctggggggcagctgctgctgcgccccgalttagctccgaattgtggagagctgggagtgctgagcagcattccagccctggcaggcggcctgtgtacatttcagcacttccagtgctggggcgcacctgtgtgcaccgcagtggtgctcacagctgctcatgtcatcagcgacaattaccagctctggctgggtcgcccaactgttttgacgacgaaaaacacagcccagttgttcatgtcagtgagagcttcccacacctgcttcaacatgagctcctctggagaaacacaccgccgaagcagcagagactacagccacgacctatgctgctccgcctgacagagctcctgtatcacatcacagagcctgtgtgaagctcgtggagctggccaccaggaaaccgaagtggggagcagctgtgttgcttccgctgtggggcagcaltcgaaccagagaattctcatlcttcagatgaltccagtgctgtggacctcaaaaltctgctctaatgtatgagtgcaaaaaagtcacgtccagaagctgacagacttcatctgtgtgtcggcacctggaggtggcgaagaacatcctgtatgtgagtgattcagctggggccctgatgtgtgtgagtgcttccaggtgtcacatcatggggctacgtccctgtgtggcccccaataagccttctgtccgcctcagagtgctccttctgtgaagtgcattcaggacgacacatgacacatgacacatgctcga

More info: Please contact us for any question: contact@proteogenix.fr