

SSX4-synovial sarcoma, X breakpoint 4 Gene

Reference: PTXBC005325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein SSX4; CT5.4; cancer/testis antigen 5.4; synovial sarcoma, X breakpoint 4; SSX family member 4		

For research use only.

GENE DESCRIPTION	synovial sarcoma, X breakpoint 4
GENE ID	6759
ORIGIN SPECIES	Human
ACCESSIONS	BC005325
VECTOR	pDONR223
NCBI SYMBOL	SSX4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggagacgacgccttgcgaaggagaccaggatgatgctcaaatcagagaaagttacgaaaggccttcgatgatattgccaatactctctaaagaagatgggaaaagatgaaatcctcggagaaaatcgtctatgtgtatatgaagctaaactatgaggtcatgactaaacta
gggttcaaggtcaccctccacccttcatgctagtaaacgggctgcagacttcacgggaatgatttggtaacgatcgaaccacagggaatcaggttgaacgtcctcagatgacttggcagcctccagagaatcttcccgaagatcatgcccaagaagccagcagaggaagaaatg
gtttgaaggaaatgccagaggcatctggcccaaaatgatgggaacagctgtgccccgggaaatccaagtaccttggagaagatcaacaagacatctggacccaaaagggggaacatgcctggaccacagactgcgtgagagaagcagctgggtgttatgaagagatca
gcgaccctgaggaagatgacgagtaa