

NUPL2-nucleoporin like 2 Gene

Reference: PTXBC005327

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CG1; NLP-1; NLP_1; hCG1; nucleoporin-like protein 2; H_RG271G13.9; NUP42 homolog; nucleoporin hCG1; nucleoporin-like protein 1; nucleoporin like 2		

For research use only.

GENE DESCRIPTION	nucleoporin like 2
GENE ID	11097
ORIGIN SPECIES	Human
ACCESSIONS	BC005327
VECTOR	pDONR223
NCBI SYMBOL	NUPL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccattgtcaattctcctcaaggccggtgctgagatcgggtctggaacgaacatccgggtctaggggtgcaggaggagggacggcagcaaccgcagcagccttcaggtaataatagcgtggatggaatacaactagccagagatattccaatgtcatcagccatccagtttccaaatccacacatggggggcagcagagatcaagaaaagccatattcagttctttgattctggagcttcaactaacaggaaggagggcttggattgtctgagaacccattgcttcacttagctctgatgagcagaagatgaaagaaactctggaaggaaattgtaaaagatatggaggtttgggaatcatcaggcagtgatgtttctgtttatcaccagtgaagaaagaaacctaattattcagggtttacagacattcaccagaggaattgaggcttgaataccataacttctaaccagcaataacttacagagttatctaaattctgccaacgtttaataatcaatggaggaacagggtaaatgaactgaaagtctaaatatatacaactaaagtagctttgaaatgccttggagatatataccacctgtctccaaagatctatggtcccgactgtgctctga