



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TFPI2-tissue factor pathway inhibitor 2 Gene

Reference: PTXBC005330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PP5; REF1; TFPI-2; placental protein 5; retinal pigment epithelium cell factor 1		
For research use only.			

GENE DESCRIPTION	tissue factor pathway inhibitor 2
GENE ID	7980
ORIGIN SPECIES	Human
ACCESSIONS	BC005330
VECTOR	pDONR223
NCBI SYMBOL	TFPI2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccacccgtctgcgccctcggggctgtcgaattctgtcgtcttctcgtcgaaggagctgcactggcgcatgtctcagcaggccaacaggaaataacgcggagatctgtctctgccctagactacggacctgcggggccctactctctcgttactactacgacaggtacacgcagagctgcgcacgaattctgtacggggctcgcggcgcaacgccacaattctacactggaggtgcgcagatcgtctgtggagtagaaaaaagtccccaaagtgtccggcgctcgaatgtagtgtgacgacacagctgtgtaggggtccacagaanaattttcttaatacgaattccatgacagtctgtgaaaaaattctttcgtgggtgtgcacccgaaggaacaggtttccagcagagctactgtatgtgctctgcgcacccaagaagaataatctcatatttctgtacgactccaaagtagggagactgtctgtcctccaaattgtactcgtcgtcgtccaaattgaatgaataatccaaatatacagaacctgtgactgtttccactactacgtctgtggaggaatgacacatttggtagcaggagatgtcaaacgtgcagtgtccaaaggctgtgaaaaaagaanaaaatgcacaaagctgcgtctccagtagggcgcggaaatttcgaagaagcaatttaa