

GTF2H2-general transcription factor IIH, polypeptide 2, 44kDa Gene

Reference: PTXBC005345

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTF2; BTF2P44; T-BTF2P44; TFIIH; p44; general transcription factor IIH subunit 2; BTF2 p44; TFIIH basal transcription factor complex p44 subunit; basic transcription factor 2 44 kDa subunit; general transcription factor IIH polypeptide 2; general transcri		

For research use only.

GENE DESCRIPTION	general transcription factor IIH, polypeptide 2, 44kDa
GENE ID	2966
ORIGIN SPECIES	Human
ACCESSIONS	BC005345
VECTOR	pDONR223
NCBI SYMBOL	GTF2H2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgaagaacctgaagaactaagcgatgggaaggaggctatgaagaacatgggagattctaaagaagatgaatctggatcactaaagctacaatagaagacattctattcaaggcaagagaaaaagagtattgagcaccatggacaagttcgacttggatgatgcgccacc
tttatgtgtagtagatggatcaagaacaatggaagaccaagatttaagcctaataagactgacgtgtactttaagttgttggaatacttttagagggaatattttgatcaaaatcctattagtcagattggaataattgtaactaagagtaaaagagctgaaaaattgactgaacttcaggaaacc
caagaaacataaacgtctttgaagaagctgtggatgacctgcatggagagccatctcttataatccctaagcatagctatgcagactctaaagttagtattatattatgtataattga