

APOC2-apolipoprotein C-II Gene

Reference: PTXBC005348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APO-CII; APOC-II; apolipoprotein C-II; apolipoprotein C2		

For research use only.

GENE DESCRIPTION	apolipoprotein C-II
GENE ID	344
ORIGIN SPECIES	Human
ACCESSIONS	BC005348
VECTOR	pDONR223
NCBI SYMBOL	APOC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcacacgactctcccagctctgtttctgtcctctgtattgggatttgaggtccaggggacccaacagcccccagcaagatgagatgcctagcccgacccctctcaccagggtgaaggaatctctccagttactgggagtcagcaagacagccgccagaacctgtacgaga
agacatacctgcccgtctgtagatgagaaactcagggaactgtacagcaaaagcacagcagccatgagcacttacacaggcatctttactgaccaagttcttctgtctgaaggagaggagtaa