

RPLP2-ribosomal protein, large, P2 Gene

Reference: PTXBC005354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D11S2243E; LP2; RPP2; 60S acidic ribosomal protein P2; acidic ribosomal phosphoprotein P2; renal carcinoma antigen NY-REN-44; ribosomal protein, large, P2; ribosomal protein lateral stalk subunit P2		

For research use only.

GENE DESCRIPTION	ribosomal protein, large, P2
GENE ID	6181
ORIGIN SPECIES	Human
ACCESSIONS	BC005354
VECTOR	pDONR223
NCBI SYMBOL	RPLP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgctacgtcgctcctacgtctggctgccctaggggcaactcctccccagcgccaaggacatcaagaagatctggacagcgtgggtatcgaggcggacgacaccggctcaacaagggtatcagtgagctgaatggaaaaacattgaagacgtcattgccagggtattgg
caagcttgcagtgatcctgctggtggggctgtagccgtctctgctgccccagcgtctgcagccctgctgctggttctgccccctgctgcagcagaggagaagaagatgagaagaaggaggagctgaagagtcagatgatgacatgggatgttgcccttttgattaa