

EGLN1-egl nine homolog 1 (C. elegans) Gene

Reference: PTXBC005369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf12; ECTY3; HALAH; HIF-PH2; HIFPH2; HPH-2; HPH2; PHD2; SM20; ZMYND6; egl nine homolog 1; HIF-prolyl hydroxylase 2; egl nine-like protein 1; hypoxia-inducible factor prolyl hydroxylase 2; prolyl hydroxylase domain-containing protein 2; zinc finger MYND		

For research use only.

GENE DESCRIPTION	egl nine homolog 1 (C. elegans)
GENE ID	54583
ORIGIN SPECIES	Human
ACCESSIONS	BC005369
VECTOR	pDONR223
NCBI SYMBOL	EGLN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttgcttgatccggcaatggaacgggttatgtacgtcatgttgataatccaaatggagatggaagatgtgtgacatgtatatattatcttaataaagactgggatgccaaggtaagtggaggtatactcgaattttccagaaggcaagcccagttgtgacattgaacccaaattgat
agactgctgtttttctgctgacgtcgcaaccctcatgaagtacaaccagcatatgctacaaggtacgcaataactgtttggtattttgagtcagatgagagagcagcagctaaagtaaaatafctaacaggtgaaaaagggtgaggggtgaactcaataaaccttcagattcggtcggtaaa
gacgtcttctag