

CAPNS2-calpain, small subunit 2 Gene

Reference: PTXBC005397

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	calpain small subunit 2; CSS2; calcium-dependent protease small subunit 2		
For research use only.			

GENE DESCRIPTION	calpain, small subunit 2
GENE ID	84290
ORIGIN SPECIES	Human
ACCESSIONS	BC005397
VECTOR	pDONR223
NCBI SYMBOL	CAPNS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttctgcaaggctctattggaaggagcagatcaggtcttgagaagctcttgaggcctcttgaggaggtggtcagagaagagaaggaggaggaagaatattggaggagatgttgagggaattgtgaattttatcagtgaggctgcagcagctcagtatactcagaaccgcc
tcccactcagcagcatttcaccagtggtgaggcctcagaaagttaggagagtttaggcgatttcggcaacaatttacacagctggctggaccagacatggaggtgggtgccactgatctgatgaattatctcaacaagtccttctaagcacaagaatcttaagactgacgggttttagcttgaca
cctgccggagcattgtgtctgcatggacagtgacacgactggttaagctggcctttgaagaatttaagtatctgtggaacaacatcaagaaatggcagtggtgttataagcagtagacaggacattctgggtctctgggaagttctcagctgcggggagctctgcaggccgagcctcc
agctaaatgaacaactttaccaatgattgtccggtatgctaatgaagatggagatatgattttaacaatttcacagctgctgtgctccgctggatgcatgtttcgtgccttcaagtcctctggatagatagagatggcctgattcaagtgctctatcaagaatggctgcagttgacctgt
attcctga