

## CENPK-centromere protein K Gene

Reference: PTXBC005400

|                |                                                                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | AF5alpha; CENP-K; FKSG14; P33; Solt; centromere protein K; SoxLZ/Sox6-binding protein Solt; interphase centromere complex protein 37; leucine zipper protein FKSG14; protein AF-5alpha |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | centromere protein K                                                                              |
| <b>GENE ID</b>          | 64105                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC005400                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CENPK                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaatcaggagatctagatccggatagtactacagatgtgggagatgttacaataactgaagaagaacttattagagaatgtgaagaatgtggaagatgtggaagaatgtcagaataaattatcatttgggaactgaaacactcaccgattcaaatgctcagctatcattgttaattatgc  
aagtaaaatgtttaaccgctgaactcagtcagtcaatggcagaaaaaacacctgaacaattcccttgactgaagacgttctcataacattaggaagaagaaggttccaaaagctgagacaagatctgaaatgggtactgtccactaaggagtcagaagaatgaaaagttaaagggaagacttagaaa  
gggaacaacgggtgggtggaacagcaacagataatggaatctttaaagtactacacagtgaaatgaaaaataaggttgaacattttctgaatcaagaatctttaaagaactgaaactaaaatgcttaataataaagaataaaggagaactcttgatgacttggcgagttctagaaga  
ccattttctctgctgatagaagtgttaaaagaaaaaagaaacattcaagaatcatctgtaaacctgataacactgcatagaatgttagagattctataaatagattttgatgtccacatgatccatgtcaaaattagtgattcctttggccacctatgttgagctgctgctgctaatgg  
aattgccttgagacatccagaagatccaaccgaataagattagaagctttccatcagtaa