



Reference: PTXBC005811

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PDHK2; PDKII; pyruvate dehydrogenase kinase, isozyme 2; PDH kinase 2; [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 2, mitochondrial; pyruvate dehydrogenase kinase, isoenzyme 2; pyruvate dehydrogenase, lipamide, kinase isozyme 2, mitochon		

GENE DESCRIPTION	pyruvate dehydrogenase kinase, isozyme 2
GENE ID	5164
ORIGIN SPECIES	Human
ACCESSIONS	BC005811
VECTOR	pDONR223
NCBI SYMBOL	PDK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1