

POLR2L-polymerase (RNA) II (DNA directed) polypeptide L, 7.6kDa Gene

Reference: PTXBC005903

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RBP10; RPABC5; RPB10; RPB10beta; RPB7.6; hRPB7.6; DNA-directed RNA polymerases I, II, and III subunit RPABC5; DNA-directed RNA polymerase III subunit L; RNA polymerase II 7.6 kDa subunit; RNA polymerases I, II, and III subunit ABC5; RPB10 homolog; polymer		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide L, 7.6kDa
GENE ID	5441
ORIGIN SPECIES	Human
ACCESSIONS	BC005903
VECTOR	pDONR223
NCBI SYMBOL	POLR2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcatccctgtacgtcttcacttgtggcaagatcgtcggcaacaagtgggagcttacctggggctgctgcaggccgagtagcaccgagggggatgcgctggatgccctgggcctgaagcgtactgctgccgccgatgctgctggccacgtggacctgacgagaagctgct
caattatgcacccctggagaagtga