

### SSX3-synovial sarcoma, X breakpoint 3 Gene

Reference: PTXBC005904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Kruppel-associated box containing gene product SSX3; protein SSX3; CT5.3; cancer/testis antigen 5.3; synovial sarcoma, X breakpoint 3; SSX family member 3		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	synovial sarcoma, X breakpoint 3
<b>GENE ID</b>	10214
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC005904
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	SSX3
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgaacggagatgacaccttgcaggagaccacggttgctcaataccagagaagatacaaaagccttcgatgatattgccaatactctctaaggaagatgggaaaagatgaaagtctgggaaaatcgtctatgtgtatatgaagaagaatgaggccatgactaaact  
aggtttcaaggccatcctccatcttcatgcftaataaacgggtcagactccaggggaatgatttgataatgaccctaaccgtgggaatcagggtcaacgtcctcagatgactttggcaggctccagggaaatctccgaagatcatgccaagaagccagaggaagaaatgt  
ttcgaaggaagtgccagaagcatctggcccacaaaacgatgggaaacagctgtgccccgggaaaaccaactacctctgagaagattaacatgatctgtgtcctcagagatactgcaggttcggtccagaccactgcaataa