



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MIA-melanoma inhibitory activity Gene

Reference: PTXBC005910

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CD-RAP; melanoma-derived growth regulatory protein; melanoma inhibitory activity

For research use only.

GENE DESCRIPTION	melanoma inhibitory activity
GENE ID	8190
ORIGIN SPECIES	Human
ACCESSIONS	BC005910
VECTOR	pDONR223
NCBI SYMBOL	MIA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgggctcctgtgtgtccttgggtgatcatcattctgtctgtcctctcgggaactgggtgtcagggggggctctatgcccaagctggctgacccggaagctgtgtgtcggacaggaggtgcagccaccctatctcatgtgtgtggcccttcaggactaatgtcccccgactgcgattc
ctgacattaccctccggccaaagtgtgtgtgtctctccaagtgcgaagccgtggggcctctctctggggcagcagcgttcaggggagattatctaggagatctgctgtcgtcgtggctatttccccagtagtcattgtccgagaggaccagaccctgaaacctggcaaaigtcatgtga
agacagacaaatggggattctctatctccagtga