

SCP2-sterol carrier protein 2 Gene

Reference: PTXBC005911

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NLTP; NSL-TP; SCP-2; SCP-CHI; SCP-X; SCPX; non-specific lipid-transfer protein; propanoyl-CoA C-acyltransferase; sterol carrier protein X		

For research use only.

GENE DESCRIPTION	sterol carrier protein 2
GENE ID	6342
ORIGIN SPECIES	Human
ACCESSIONS	BC005911
VECTOR	pDONR223
NCBI SYMBOL	SCP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtttccggaagccgccagttctttagaactcatcaattgaagctgttccaaccagctctgcaagtgatgatttaaggcaaatctgttttaaggagattgagaagaactgaagaggaaggggaacagttgtgaagaaatcggtggtattttgccttcaagtgaaagatggcc
ctgggggtaagaggccacctgggtggtgatgtgaagaatggcaaggatcagtgcttctaactcagataagaaggctgactgcacaatcaaatggctgactcagacttctgctttaatgactggtataatgaatcctcagtcggccttcttcaaggcaaatgaaatcactggca
acatgggtctcgtatgaagttacaaaatcttcagcttcagccaggaacgctaagctctga