

PLA2G2A-phospholipase A2, group IIA (platelets, synovial fluid) Gene

Reference: PTXBC005919

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOM1; PLA2; PLA2B; PLA2L; PLA2S; PLAS1; sPLA2; phospholipase A2, membrane associated; GIIC sPLA2; NPS-PLA2; group IIA phospholipase A2; non-pancreatic secretory phospholipase A2; phosphatidylcholine 2-acylhydrolase 2A; phospholipase A2, group IIA (platele		

For research use only.

GENE DESCRIPTION	phospholipase A2, group IIA (platelets, synovial fluid)
GENE ID	5320
ORIGIN SPECIES	Human
ACCESSIONS	BC005919
VECTOR	pDONR223
NCBI SYMBOL	PLA2G2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagacctctactgttggcagtgatcatgatttggcctactgcaggcccatgggaatttggtaattccacagaatgatcaagttgacgacaggaaggagccgcactcagttatggcttctacggctccactgtggcgtgggtggcagaggatcccccaaggatgcaacgga
tcgctgctgtgtcactcatgactgttctacaacgtctgggagaacgtggatgtggcaccaatttctgagctacaagtttagcaactcggggagcagaatcacctgtgcaaacaggactcctgcagaagtcactgtgtgagtgatgataaggctgtgccacctgttttctagaaacaag
acgacctacaataaaaagtaccagtactattccaataaacactgcagaggagcaccctcgttgcctga