

NIPSNAP3A-nipsnap homolog 3A (C. elegans) Gene

Reference: PTXBC005935

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC299; NIPSNAP4; TASSC; protein NipSnap homolog 3A; protein NipSnap homolog 4; target for Salmonella secreted protein C; nipsnap homolog 3A		

For research use only.

GENE DESCRIPTION	nipsnap homolog 3A (C. elegans)
GENE ID	25934
ORIGIN SPECIES	Human
ACCESSIONS	BC005935
VECTOR	pDONR223
NCBI SYMBOL	NIPSNAP3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtcctcagaagcgccctgactcggcgctggcctcagcgagctggcgctcagatgtgctcatcttttctacgggaccagacaatacagtggaattatctatgaatttcgtcttattacctaagccctcaaatgaatgagttcctggaaaatttgagaaaaacgctcatctc
ggacagctcactctgaattggttgatactggagtgtagaatttggaggcagaatgaatacaggttttcatatttggagtgatgataatttgcctcatcgaactgaagttcggaaagccttggccaagataagggaatggcaagaacaattctcattccaaatttggctcattgataacaaga
gagtgagattactatctgtgacatggtgcaaatagaaaaacccaaaagaaggagtgatgaactggccacttttcagatgaacctgtgtggccagctctgtgggtgatgcatttaaaaggcagttcatgctcatgtaactaggtacacaaaactagttggagttccacacag
agtacggagcactcaacagagttcatgtcttttggggaatgagagtgacagatagtcgtgcagctgggagacataagtcctcagagatccagagttgtggcagctgttcgggaaagtgtaactacctaagatctcagcagaatatgctctgattcctacatcgtttcaccactgaaatag