

NDUFA12-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 12 Gene

Reference: PTXBC005936

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B17.2; DAP13; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 12; 13 kDa differentiation-associated protein; CI-B17.2; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 12; NADH-ubiquinone oxidoreductase subunit B17.2; complex I B17.2 subunit		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 12
GENE ID	55967
ORIGIN SPECIES	Human
ACCESSIONS	BC005936
VECTOR	pDONR223
NCBI SYMBOL	NDUFA12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtttagtcaggctcctgaacgcgggctgcagcagatcaccggccacggcgggtctccgaggtatctacgggttttttcaggacaaatgatgcgaaggttggtacattagtggtgggaagacaaatgatgaaacaaatactatgaagacaacaagcaatttttgccgtcaccgatg
gggtgtatatactactgaaatgaatggcaaaaacacattctgggatgtggatggaagcatgggtgcctctgaatggcatcgttggcttcacagtafagctgatgatcctccaacaacaaaccattgctgctcgttaattcatttggacgaaccataaattcaacgtgactggcaccacagaac
aatatgtaccttattctaccactagaagaagattcaggagtggtatccacctcaacaccttacaagtaa