

SMPX-small muscle protein, X-linked Gene

Reference: PTXBC005948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DFN6; DFNX4; small muscular protein; deafness, X-linked 6, sensorineural; stretch-responsive skeletal muscle protein; small muscle protein, X-linked		

For research use only.

GENE DESCRIPTION	small muscle protein, X-linked
GENE ID	23676
ORIGIN SPECIES	Human
ACCESSIONS	BC005948
VECTOR	pDONR223
NCBI SYMBOL	SMPX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatatgtcgaacagccagtttccaatgttagagccatcaggcaaatatcaatattccaatgggagccttccggccaggagcaggtcaacccccagaagaaaagaatgtactctgaagtggaggagggtgttctccacctcgatgaggagaagaagccaattccaggagcg
aagaaacttcaggacctgcagtcattctatcggaatccagaatattaaaagtgaaactaaaatattgtcccaaaagctgaacagtag