

RLN1-relaxin 1 Gene

Reference: PTXBC005956

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1RLX; RLXH1; bA12D24.3.1; bA12D24.3.2; prorelaxin H1; preprorelaxin H1; relaxin 1		

For research use only.

GENE DESCRIPTION	relaxin 1
GENE ID	6013
ORIGIN SPECIES	Human
ACCESSIONS	BC005956
VECTOR	pDONR223
NCBI SYMBOL	RLN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctcgcctgttctgtccacctgctagaattctgtttactactgaaccaattttccagagcagtcgcccgaatgaagacgatgttataaattatcgccgcgaattagtcgcgcagattgccatttgcggcatgagcacctggagcaaaaggctctgagccaggaagatgct
cctcagacacctagaccagtggcagaaattgtaccatcctcatcaacaagatacagaaactataattatcatgttggaattcattgctaatttgcaccggagctgaagcgagccctatctgagaggcaaccatcattaccagagctacagcagtatgtacctgcaataaggattccaatctt
agcttgaagaatttaagaacttattcgaataggcaagtgagccgcagacagcaatccttcagaattaaatacttaggcttggatactcattctcaaaaaagagacgacctacgtggcactgttgagaatgttgcccaattgggtgtacaaaaggctctctgctaaatattgctga