

RHOA-ras homolog gene family, member A Gene

Reference: PTXBC005976

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	small GTP binding protein RhoA; transforming protein RhoA; ARH12; ARHA; RHO12; RHOH12; Aplysia ras-related homolog 12; oncogene RHO H12; ras homolog family member A		

For research use only.

GENE DESCRIPTION	ras homolog gene family, member A
GENE ID	387
ORIGIN SPECIES	Human
ACCESSIONS	BC005976
VECTOR	pDONR223
NCBI SYMBOL	RHOA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccatccggaagaactggtgattgttgatggagcctgtgaaagacatgcttctcatagcttcagcaaggaccagtcccagaggtgatgtgccacagtgtttgagaactatgtgcagatcagagtgatggaagcaggtagaattggctttgtggacacagct
ggcgaggaaagattatgatcctgagggccctctctaccagataccgatgttatactgatgtgttttccatcgacagcctgatagtttagaaaacatcccagaaaagtgaccaccagaagtcaagcatttctgtcccaactgcccatacctgttggataagaaggatctcgaa
tcatgagcacacaaggcgggagctagccaagatgaagcaggagccggtgaaacctgaagaagcagagatatggcaaacaggattggcctttgggtacatggaggttcagcaagaccaagatggagtgagagaggttttgaattggctacgagagctctgcaagctaga
cgtgggaagaaaaatctgggtccttctgtga