

PHF6-PHD finger protein 6 Gene

Reference: PTXBC005994

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BFLS; BORJ; CENP-31; PHD finger protein 6; PHD-like zinc finger protein; centromere protein 31		

For research use only.

GENE DESCRIPTION	PHD finger protein 6
GENE ID	84295
ORIGIN SPECIES	Human
ACCESSIONS	BC005994
VECTOR	pDONR223
NCBI SYMBOL	PHF6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgtcaagctcagttgaacagaaaaaggcctacaagacagcgcgaatgtggcttttgaagcaaatagagacaagaatgtggacagttactaatctgaaaccagaagtgccagcgcaccataagtcagctctttcatctgctttggtatcatcacacttgataatgaaagtctt
ggfaggattttattgaagatgtccaaaaggaaftaaaagaggcagcaagctgatgttctttgfccttctctggagcaacaattggtgtgatgtgaaacatgtcacaggacataccactaccactgtgcatgcatgataaagctcaaatcagagagaaacctcacaagaaattac
atggcctatgccgaaaacacaagaaaactgcacataactcgaagcagctgattgaagaagaatttaataacatgaactggagccctcaccctaaaagtaaaaagaaaagtcgaaaggaagccaaagaaaactaatttaagggtgtcagaagataccaggtccacatcctc
ccatggaacagatgaaatggaagtagttctatagataggtctccacagaaagcgccttagtgacaccaggcctaattgtgallttgccatgtaggggaggaagaaatgaagcacgaggaactgcataatlaagccaagaagcagctgccattataagtcagttgtt
tttctctggcagctccagctcacaacaacatcaagagcagaalltggagactttgatataaaactgtacttcaggagattaaacgaggaagaaagaaatggctctgtagttttatattgttatgcaacattacactgatttgcctttaaatttagagataccctcaaatatccagtcacagaa
attaaagtag
```