

C6orf125-chromosome 6 open reading frame 125 Gene

Reference: PTXBC006007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf125; Cbp6; M19; MNF1; bA6B20.2; ubiquinol-cytochrome-c reductase complex assembly factor 2; breast cancer-associated protein SGA-81M; cytochrome B protein synthesis 6 homolog; mitochondrial nucleoid factor 1; mitochondrial protein M19; ubiquinol-cyto		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 125
GENE ID	84300
ORIGIN SPECIES	Human
ACCESSIONS	BC006007
VECTOR	pDONR223
NCBI SYMBOL	C6orf125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggccagccggtaccggcgtttcttaagctctgtgaggaatggccagtggacgagaccaaacggggccgggacttggcgcttacctgcgacagcgggtagcacaggccttccgggaggagagaataaccagggttcagagcctgaggcctgtgatcagatgtacgagagc
ttagcgcgactccattcaaaactactacaacacaagtacctcgccccagagacaccagcttcagtgccctgtcgttggagagtacaagctgatcctgtccacagacaccttggagagccttaaggaaatagataaaggcatgtggagaactgcaggagaagtttgcaccaagggtc
ctgaggagatcataaggcctga