

DDI2-DDI1, DNA-damage inducible 1, homolog 2 (*S. cerevisiae*) Gene

Reference: PTXBC006011

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNA-damage inducible protein 2 (DDI2); protein DDI1 homolog 2; DDI1, DNA-damage inducible 1, homolog 2; DNA damage inducible 1 homolog 2		

For research use only.

GENE DESCRIPTION	DDI1, DNA-damage inducible 1, homolog 2 (<i>S. cerevisiae</i>)
GENE ID	84301
ORIGIN SPECIES	Human
ACCESSIONS	BC006011
VECTOR	pDONR223
NCBI SYMBOL	DDI2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtcaccgtgtactgtgtgctggaggacctctccgaggtgacctttccctccaggtcgacgacgacttcgagctgcacaaactccgcgctgtgctgagctgagctgtggcatccccgcagccgagagccagatctctatgctggaagacctctcacagacaaccacagatcat
tggcttcttatgcttgaagatggggacgtgtgtattttacgacagaaggagaaatgcagacctcgacctccagtgacgttcccaacttaccggaatagatttcagtagtatagctgtgctggcacatcaagtccccggcagcgccagccaggaacacagcagtcctcactcatctc
ctggagaataaactcatctctcagggcttggacaatccagccttgcctccgagataatgtgtgctggccaacccgcatgagctgtcctgtctgaaggaaacgcaatccacccctggcagaagctctgctcagtgaggaccttgagaaattttctagagctctgtgtggagcagcagcaggaccga
gccccggagagagcaagaaggattcgtctgtttctgtgatcccttgaccttgaagctcaggcagaagatagaagaatataaaggttaa