



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ORAI3-ORAI calcium release-activated calcium modulator 3 Gene

Reference: PTXBC006126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM142C; protein orai-3; transmembrane protein 142C; ORAI calcium release-activated calcium modulator 3		

For research use only.

GENE DESCRIPTION	ORAI calcium release-activated calcium modulator 3
GENE ID	93129
ORIGIN SPECIES	Human
ACCESSIONS	BC006126
VECTOR	pDONR223
NCBI SYMBOL	ORAI3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]