

C19orf48-chromosome 19 open reading frame 48 Gene

Reference: PTXBC006151

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C19orf48; multidrug resistance-related protein; chromosome 19 open reading frame 48		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 48
GENE ID	84798
ORIGIN SPECIES	Human
ACCESSIONS	BC006151
VECTOR	pDONR223
NCBI SYMBOL	C19orf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagtgtggaggccgtcttgagatccaggccatcactggcagcaggtgtctctcatgtgtccaggcccgccaggccaccaggtcatgtgggaccaaccagtgacaaaggacttggtgctgagccacacccaggagaaggtggataagtggtaccagggttctctgaggctaggggaggagccaccccgcttcctattgtgaccaggccatggggaggagctgtccatacgccaccgtgagacctgggcctggctctcaaggacagacaccgcctggcctggtgctccagggtgaagcaggccagaatcctggggagctgctcctggtttga