



Reference: PTXBC006155

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADTAB; CLAPA2; HIP-9; HIP9; HYPJ; AP-2 complex subunit alpha-2; 100 kDa coated vesicle protein C; adapter-related protein complex 2 subunit alpha-2; adaptin, alpha B; adaptor-related protein complex 2 subunit alpha-2; alpha-adaptin C; Huntingtin interacti		

GENE DESCRIPTION	adaptor-related protein complex 2, alpha 2 subunit
GENE ID	161
ORIGIN SPECIES	Human
ACCESSIONS	BC006155
VECTOR	pDONR223
NCBI SYMBOL	AP2A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]