



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DCTN1-dynactin 1 (p150, glued homolog, Drosophila) Gene

Reference: PTXBC006163

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DAP-150; DP-150; P135; dynactin subunit 1; 150 kDa dynein-associated polypeptide; dynactin 1 (p150, glued homolog, Drosophila)		

For research use only.

GENE DESCRIPTION	dynactin 1 (p150, glued homolog, Drosophila)
GENE ID	1639
ORIGIN SPECIES	Human
ACCESSIONS	BC006163
VECTOR	pDONR223
NCBI SYMBOL	DCTN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacggccggccagctcgctggaagagatccacacgtctgtctcagcagatctctgcgatgggctgcacatctccagctccagcatgagagaacacatctccaaggagcccgatgaagcgcactctggcatccctccccctctcactgtttgcaaaagctatcccatgaggccctggc
agtgaggtaccagctcgaggcgctgtatctgaagacagccagctcgtctggagacattgaatcaattgagcacacacgcacgtagttagatcacatcctgcacaggccctgctgccaaagagccctgcgcgcccaacttaaggagcaatggctcagcttaagtccttgagtgacacctgcgag
agctcaacagatgaggtctccaaggagacagatgatctgcgcctcgaggccagatgaccacatgactttgccaccttccctctcatcagctctctcagggccaagagagcagcatgagatgacacagctcatcagggccaaagtgacatctctcatgctcgctgtttggacagcgcacacc
ggctctctgcacccggacgcagctccacacgtctcacagtctcagatctcctcaaa