

ZNF397-zinc finger protein 397 Gene

Reference: PTXBC006172

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF47; ZSCAN15; zinc finger protein 397; zinc finger and SCAN domain-containing protein 15; zinc finger protein 47		

For research use only.

GENE DESCRIPTION	zinc finger protein 397
GENE ID	84307
ORIGIN SPECIES	Human
ACCESSIONS	BC006172
VECTOR	pDONR223
NCBI SYMBOL	ZNF397
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcttggaatctggagtgattcaaccctgatacctcagatcctccggaacaagaactaactagtgaagtagaagataactttctgggatgagaatttaagcagaatgggagctactcaatcctccaagaattgttcgtcagaattcagaataattgctaccaggagacact
gggccccgggagctctgagccgactccaggaaacttgcatacgtgctaatgccagagttgcacacaaaggagcagatctagaactcgtgctactggagcagttctgagcattctgctgaggagctgcagatctgggtcagcaacataatccagaagcggcgaggaagctfg
acctgttgaggattagagagggagttgatgaccagggcagcaggtcccagctagtccacagggaccagcagtgccatggaaggatttaacatgtctcagagcatccaagagtcacaagacatccacctccagccctaaagacacagctgaaatctggaacctgccttcc
cctaaaagtattgtgagaacagtgaaacagcaacaaaaggggcatcagaagaaaaatcacagggactcctcaggaaacctcatttcagggaattaagttgtccagacctccaaagctctcagctatccgtgggaatgtttctcaggaagtttccggcgatcatagctgc
tgaggctacacatcaacaattctgtttgcatcaacacttgcagccacctcctgcatctaaaaatgtaataagaagatattttctga
```