

INF2-inverted formin, FH2 and WH2 domain containing Gene

Reference: PTXBC006173

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf151; C14orf173; CMTDIE; FSGS5; pp9484; inverted formin-2; HBEAG-binding protein 2 binding protein C; HBEBP2-binding protein C; inverted formin, FH2 and WH2 domain containing		

For research use only.

GENE DESCRIPTION	inverted formin, FH2 and WH2 domain containing
GENE ID	64423
ORIGIN SPECIES	Human
ACCESSIONS	BC006173
VECTOR	pDONR223
NCBI SYMBOL	INF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggtagaaggagggcgccagcgcgaagtgggcagcgtgaaggagaagctggggccacaggattcggaacccacggaggccaacctggagagcgcggacccctgagctgtgcatccggctgctccagatgccctctgtgtcaactactcggcctgcgaagcgctggag
ggcagcgcagggcgctggatggatgagttcctggagcagagcggcctggacctgctgctggagggcgtggcgcggtgtcgggcccggcggttgcaactatctccgacgacctgctgcagctcacctgcgtcagctgctgctggcgccgctcatgaactcggcgaggcagtcagta
catcctcagcaaccagggtctgctgcccagctctccaggccctggacacatcaacgtgatggtagaagaagcagggtgttgagctactggtgcccctgtgcatctactctcccaggggccacgtgctgacctggagcgcctggaccactacaagcgggtgtgagccagcagta
gcttcagcattgtcatgaacgagctctccggcagcgacaacgtgccctacgtgggtcacctgcttagcgtgatcaacgccgtcatcttgggccccgaggacctgcgcgcgcgacccagctgcggaacgagtttatcggtgagctgctggacgtctggctcgctcggtga