



Reference: PTXBC006175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Jun proto-oncogene, AP-1 transcription factor subunit; v-jun sarcoma virus 17 oncogene homolog; v-jun avian sarcoma virus 17 oncogene homolog; proto-oncogene c-Jun; jun proto-oncogene; jun oncogene; Jun activation domain binding protein; AP-1; AP1; transc		

For research use only.

GENE DESCRIPTION	jun oncogene
GENE ID	3725
ORIGIN SPECIES	Human
ACCESSIONS	BC006175
VECTOR	pDONR223
NCBI SYMBOL	JUN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1