

TNFRSF9-tumor necrosis factor receptor superfamily, member 9 Gene

Reference: PTXBC006196

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4-1BB; CD137; CDw137; ILA; tumor necrosis factor receptor superfamily member 9; 4-1BB ligand receptor; CD137 antigen; T cell antigen ILA; T-cell antigen 4-1BB homolog; homolog of mouse 4-1BB; induced by lymphocyte activation (ILA); interleukin-activated r		

For research use only.

GENE DESCRIPTION	tumor necrosis factor receptor superfamily, member 9
GENE ID	3604
ORIGIN SPECIES	Human
ACCESSIONS	BC006196
VECTOR	pDONR223
NCBI SYMBOL	TNFRSF9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacagctgttacaacatagtagccactctgttctgtctcaactttgagaggacaaagatcattgcaggatcctttagtaactgccagctgtacattctgtataataacaggaatcagattgcagtcctgtctccaatagtttccagcgcaggtggacaaaggacctgt
gacatagcaggcagtgtaaaaggtgtttcaggaccaggaaggagtggtctccaccagcaatgcagagtgtgactgcactccagggtttcactgcctgggggcaggatgcagcatgtgtgaacaggattgtaaacaaggtcaagaactgacaaaaaaggtgtgaaagactgttgccttgg
gacattaacgatcagaacgtggcatctgtcgaacctggacaaactgttcttggatggaaagtctgtcttgaatgggacgaaggagaggacgtgtgtgtgacctctccagccgacctctctccgggagcatcctgtgacccccgtgccctgcgagagagccaggacac
tctccgagatcatctcttcttcttgcgtgacgtcgtgtgtctctctctcctcagctccgttctctgtttaaaggggcagaaagaaactcctgtatatattcaacaaccatttatgagaccagtacaactactcaagaggagatggctgtagctgccgattccag
aagaagaagaggaggtgaactgtga