



Reference: PTXBC006212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leucine zipper putative tumor suppressor 2; leucine zipper, putative tumor suppressor 2; protein LAPSER1; leucine zipper tumor suppressor 2		

For research use only.

GENE DESCRIPTION	leucine zipper, putative tumor suppressor 2
GENE ID	84445
ORIGIN SPECIES	Human
ACCESSIONS	BC006212
VECTOR	pDONR223
NCBI SYMBOL	LZTS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]