



Reference: PTXBC006213

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	matrix-remodeling-associated protein 8; limitrin; matrix-remodelling associated 8; matrix remodeling associated 8		

For research use only.

GENE DESCRIPTION	matrix-remodelling associated 8
GENE ID	54587
ORIGIN SPECIES	Human
ACCESSIONS	BC006213
VECTOR	pDONR223
NCBI SYMBOL	MXRA8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]